



LKT Laboratories, Inc
 2233 University Avenue West
 St. Paul, Minnesota 55114 USA
 Ph: 651.644.8424 Fax: 651.644.8357
 www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	
Product ID	A0099
Product Name	A-779
Chemical Name	
Synonym	
Formula	C ₃₉ H ₆₀ N ₁₂ O ₁₁
Formula Wt.	872.99
Melting Point	
Boiling Point	
Purity	
Solubility	
Stability	
Store Temp	
Ship Temp	
Product Description	A potent and selective antagonist for the heptapeptide angiotensin-(1-7)[ANG-(1-7)]. It prevents ANG-(1-7)-induced inhibition of angiogenesis.
References	Santos, R.A; Campagnole-Santos, M.J.; Baracho, N.C., et. al. Brain Res. Bull 35:293-298 (1994). Machado, R.D.; Santos, R.A., Andrade, S.P. Amer. J Physio-Reg. Integ. & Comp. Physio 280:R994-RR1000 (2001).

H-Asp-Arg-Val-Tyr-Ile-His-D-Ala-OH

Pricing

1 mg	\$64.50
2 mg	\$109.30
5 mg	\$193.50

Bulk quantities available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
 St. Paul, MN 55114
 USA

Phone: 888-LKT-LABS (558-5227)
 651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.