




LKT Laboratories, Inc
2233 University Avenue West
St. Paul, Minnesota 55114 USA
Ph: 651.644.8424 Fax: 651.644.8357
www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	76455-82-8
Product ID	A0101
Product Name	A23187, 4 Bromo
Chemical Name	
Synonym	4-Bromo-A-23187, 4-BrA23187, 4-Bra23187
Formula	C ₂₉ H ₃₆ BrN ₃ O ₆
Formula Wt.	602.52
Melting Point	
Boiling Point	
Purity	
Solubility	
Stability	
Store Temp	
Ship Temp	
Product Description	A halogenated analogue of A23187. It is a nonfluorescent Ca ²⁺ ionophore used in calibration for determining cytoplasmic Ca ²⁺ by fluorescent probes.
References	Deber CM., Tom-Kun J., Mack E., Grinstein S. Anal Biochem 146: 349-352 (1985)



Pricing

1 mg	\$80.60
5 mg	\$329.30

Bulk quantities available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
St. Paul, MN 55114
USA

Phone: 888-LKT-LABS (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.