




LKT Laboratories, Inc
2233 University Avenue West
St. Paul, Minnesota 55114 USA
Ph: 651.644.8424 Fax: 651.644.8357
www.lktlabs.com info@lktlabs.com

Product Information

| | |
|----------------------|--|
| CAS No. | 50-78-2 |
| Product ID | A0819 |
| Product Name | Acetylsalicylic Acid / Aspirin |
| Chemical Name | 2-(Acetyloxy)benzoic acid |
| Synonym | Salicylic acid acetate |
| Formula | C ₉ H ₈ O ₄ |
| Formula Wt. | 180.16 |
| Melting Point | 136°C-140°C |
| Boiling Point | |
| Purity | ≥98% |
| Solubility | Soluble in ethanol, DMSO |
| Stability | 2 yrs |
| Store Temp | Ambient |
| Ship Temp | Ambient |

Product Description White granular powder. A non-steroidal anti-inflammatory agent. Has undergone clinical trials as a chemopreventive. Induces apoptosis.



Pricing

| | |
|-------|---------|
| 500 g | \$50.40 |
| 1 kg | \$80.64 |
| | |
| | |
| | |
| | |

Bulk quantities available upon request

References Kelloff, *et al. Cancer Epidemiol. Biomarkers Prev.* **3**:85-98 (1994).
Qiao, L., Hanif, R., Sphicas, E., Shiff, S.J., Rigas, B. *Biochem. Pharmacol.* **55**:53-64 (1998).
Castano, *et al. Br. J. Cancer.* **81**:294-299 (1999).

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
St. Paul, MN 55114
USA

Phone: 888-LKT-LABS (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com