




**LKT Laboratories, Inc**  
2233 University Avenue West  
St. Paul, Minnesota 55114 USA  
Ph: 651.644.8424 Fax: 651.644.8357  
www.lktlabs.com info@lktlabs.com

# Product Information

<b>CAS No.</b>	51037-30-0
<b>Product ID</b>	A0833
<b>Product Name</b>	Acipimox
<b>Chemical Name</b>	5-Methylpyrazinecarboxylic acid 4-oxide
<b>Synonym</b>	Neotigason; Soriatane
<b>Formula</b>	C <sub>6</sub> H <sub>6</sub> N <sub>2</sub> O <sub>3</sub>
<b>Formula Wt.</b>	154.12
<b>Melting Point</b>	177°C-180°C
<b>Boiling Point</b>	
<b>Purity</b>	≥98%
<b>Solubility</b>	Water 20 mg/mL
<b>Stability</b>	
<b>Store Temp</b>	Ambient
<b>Ship Temp</b>	
<b>Product Description</b>	Yellowish white crystalline powder. Acipimox inhibits lipolysis in peripheral tissues and induces large reduction in circulating serum-free fatty acids. The mechanism of triglyceride lowering appears to be an increase of VLDL-Tg clearance.
<b>References</b>	Nuutinen, J. Minn, H., Bergmand, J. <i>et al Br. J. Cancer</i> <b>80</b> :513-518 (1999). Hannah, J.S., Bodkin, NL., Paidi., MS., <i>et al. Acta Diabetol.</i> <b>32</b> :279-283 (1995).



## Pricing

1 g	\$68.40
5 g	\$265.80
10 g	\$455.40

**Bulk quantities available upon request**

## LKT Laboratories, Inc.

Address: 2233 University Ave. W.  
St. Paul, MN 55114  
USA

Phone: 888-LKT-LABS (558-5227)  
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

**Caution:** This product is intended for laboratory and research use only. It is not for human or drug use.