




LKT Laboratories, Inc
2233 University Avenue West
St. Paul, Minnesota 55114 USA
Ph: 651.644.8424 Fax: 651.644.8357
www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	616-91-1
Product ID	A0918
Product Name	N-Acetyl-L-Cysteine
Chemical Name	L- α -Acetamido- β -mercaptopropionic acid
Synonym	Acetylcysteine; Brunac; Fabrol; Fluimucil; Fluprowit; Mucosil
Formula	C ₅ H ₉ O ₃ S
Formula Wt.	163.19
Melting Point	106°C-108°C
Boiling Point	
Purity	≥98%
Solubility	Soluble in water.
Stability	
Store Temp	Ambient
Ship Temp	
Product Description	White powder. An antioxidant found to be effective in the early stages of carcinogenesis.
References	Bongers, V., deJong, J., Steen, I., De Vries, N., Bast, A., Snow, G.B., Braakhuis, B.J. <i>Eur J. Cancer.</i> 31A :921-923 (1995).



Pricing

10 g	\$23.30
25 g	\$40.30

Bulk quantities available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
St. Paul, MN 55114
USA

Phone: 888-LKT-LABS (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.