



**LKT Laboratories, Inc**  
 2233 University Avenue West  
 St. Paul, Minnesota 55114 USA  
 Ph: 651.644.8424 Fax: 651.644.8357  
 www.lktlabs.com info@lktlabs.com

# Product Information

<b>CAS No.</b>	
<b>Product ID</b>	A0963
<b>Product Name</b>	ACTH (1-10), human
<b>Chemical Name</b>	
<b>Synonym</b>	
<b>Formula</b>	C <sub>59</sub> H <sub>78</sub> N <sub>16</sub> O <sub>16</sub> S <sub>1</sub>
<b>Formula Wt.</b>	1299.4
<b>Melting Point</b>	
<b>Boiling Point</b>	
<b>Purity</b>	
<b>Solubility</b>	
<b>Stability</b>	
<b>Store Temp</b>	
<b>Ship Temp</b>	
<b>Product Description</b>	Adrenocorticotrophic hormone.
<b>References</b>	

H-Ser-Tyr-Ser-Met-Glu-His-Phe-Arg-Trp-Gly-OH

### Pricing

1 mg	\$28.70
2 mg	\$48.40
5 mg	\$86.00

**Bulk quantities available upon request**

### LKT Laboratories, Inc.

Address: 2233 University Ave. W.  
 St. Paul, MN 55114  
 USA

Phone: 888-LKT-LABS (558-5227)  
 651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

**Caution:** This product is intended for laboratory and research use only. It is not for human or drug use.