




LKT Laboratories, Inc
2233 University Avenue West
St. Paul, Minnesota 55114 USA
Ph: 651.644.8424 Fax: 651.644.8357
www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	59277-89-3
Product ID	A1096
Product Name	Acycloguanosine/Acyclovir
Chemical Name	Acyclovir; Vipral; Zovirax
Synonym	
Formula	C ₈ H ₁₁ N ₅ O ₃
Formula Wt.	225.20
Melting Point	256.5-257°C
Boiling Point	
Purity	
Solubility	Poor solubility in water.
Stability	
Store Temp	Ambient
Ship Temp	
Product Description	An orally active acyclic nucleoside with inhibitory activity towards several herpes viruses.
References	



Pricing

50 mg	\$61.30
100 mg	\$110.30
500 mg	\$433.20

Bulk quantities available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
St. Paul, MN 55114
USA

Phone: 888-LKT-LABS (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.