




LKT Laboratories, Inc
2233 University Avenue West
St. Paul, Minnesota 55114 USA
Ph: 651.644.8424 Fax: 651.644.8357
www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	121268-17-5
Product ID	A4515
Product Name	Alendronate Monosodium Trihydrate
Chemical Name	4-Amino-1-hydroxybutylidene-1,1-bisphosphate
Synonym	Alendronic acid; ABDP
Formula	$C_4H_{12}NNaO_7P_2 \cdot 3 H_2O$
Formula Wt.	325.08
Melting Point	233°C -235°C (dec)
Boiling Point	
Purity	≥98%
Solubility	Insoluble in water.
Stability	
Store Temp	Ambient
Ship Temp	
Product Description	A second generation bisphosphonate used as bone resorptive inhibitor.
References	O'Neil, M. (ed.) <i>The Merck Index, 13th Ed.</i> Merck & Co, Inc. Whitehouse Station, NJ 2001.



Pricing

100 mg	\$51.10
500 mg	\$178.20

Bulk quantities available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
St. Paul, MN 55114
USA

Phone: 888-LKT-LABS (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.