



Product Information

CAS Registry No. 1134-47-0

Product Name: Baclofen

Catalog No. B0110

Chemical Name: β -(Aminomethyl)-4-chlorobenzenepropanoic acid

Synonym(s): Baclon; Clofen; Lioresal

Mol. formula: $C_{10}H_{12}ClNO_2$

Formula wt: 213.66

M.p. 206°C-208°C

B.p.

Purity: $\geq 98\%$


Solubility: Soluble in water or NaOH.

Stability (in original shipping container):

Storage temp: Ambient

Product Description:

Baclofen is a GABA-b receptor agonist.



References:

Brogden, R.N., Speight, T.M., Avery, G.S. *Drugs*. **8**:1-14 (1974).

PRICING

1 g	\$26.30
5 g	\$109.80
10 g	\$175.80

*Bulk quantity available
upon request*

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.