




LKT Laboratories, Inc
2233 University Avenue West
St. Paul, Minnesota 55114 USA
Ph: 651.644.8424 Fax: 651.644.8357
www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	1134-47-0
Product ID	B0110
Product Name	Baclofen
Chemical Name	β -(Aminomethyl)-4-chlorobenzenepropanoic acid
Synonym	Baclon; Clofen; Lioresal
Formula	C ₁₀ H ₁₂ ClNO ₂
Formula Wt.	213.66
Melting Point	206°C-208°C
Boiling Point	
Purity	≥98%
Solubility	Soluble in water or NaOH.
Stability	
Store Temp	Ambient
Ship Temp	
Product Description	Baclofen is a GABA-b receptor agonist.
References	Brogden, R.N., Speight, T.M., Avery, G.S. <i>Drugs</i> . 8:1-14 (1974).



Pricing

1 g	\$29.50
5 g	\$123.00
10 g	\$196.90

Bulk quantities available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
St. Paul, MN 55114
USA

Phone: 888-LKT-LABS (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.