




LKT Laboratories, Inc
2233 University Avenue West
St. Paul, Minnesota 55114 USA
Ph: 651.644.8424 Fax: 651.644.8357
www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	132-69-4
Product ID	B1640
Product Name	Benzydamine HCl
Chemical Name	N,N-Dimethyl-3-[[1-(phenylmethyl)-1H-indazol-3-yl]oxy]-1-propanamine hydrochloride
Synonym	Benzindamine hydrochloride
Formula	C ₁₉ H ₂₃ N ₃ O•HCl
Formula Wt.	345.87
Melting Point	156°C-161°C
Boiling Point	
Purity	≥98%
Solubility	Soluble in water, ethanol,
Stability	
Store Temp	Ambient
Ship Temp	Ambient

Product Description Benzydamine HCl is a nonsteroidal analgesic and anti-inflammatory compound.

References Kim, J.H.; Chu, F.C.; Lakshmi, V.; Houde, R. **1986**. Benzydamine HCl, a new agent for the treatment of radiation mucositis of the oropharynx. *Am. J. Clin. Oncol.* 9:132-134.



Pricing

5 g	\$28.00
25 g	\$106.40

Bulk quantities available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
St. Paul, MN 55114
USA

Phone: 888-LKT-LABS (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.