




LKT Laboratories, Inc
2233 University Avenue West
St. Paul, Minnesota 55114 USA
Ph: 651.644.8424 Fax: 651.644.8357
www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	3012-37-1
Product ID	B1656
Product Name	Benzyl thiocyanate, 97%
Chemical Name	
Synonym	Benzyl rhodanide
Formula	C ₈ H ₇ NS
Formula Wt.	149.21
Melting Point	39°C-41°C
Boiling Point	230°C-235°C
Purity	≥97%
Solubility	
Stability	
Store Temp	Ambient
Ship Temp	Ambient
Product Description	Benzyl thiocyanate is effective against methylazoxymethanol acetate-induced intestinal carcinogenesis. Pale yellow liquid.
References	Sugie, S.; Okamoto, K.; Okumura, A.; Tanaka, T.; Mori, H. <i>Carcinogenesis</i> . 8 :1555-1560 (1994).



Pricing

50 g	\$22.40
100 g	\$55.40

Bulk quantities available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
St. Paul, MN 55114
USA

Phone: 888-LKT-LABS (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.