




LKT Laboratories, Inc
2233 University Avenue West
St. Paul, Minnesota 55114 USA
Ph: 651.644.8424 Fax: 651.644.8357
www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	106-51-4
Product ID	B1853
Product Name	1,4-Benzoquinone
Chemical Name	2,5-Cyclohexadiene-1,4-dione
Synonym	Quinone; <i>p</i> -quinone; 1,4-Cyclohexadienedione
Formula	C ₆ H ₄ O ₂
Formula Wt.	108.09
Melting Point	115.7°C
Boiling Point	
Purity	≥98%
Solubility	Slightly soluble in water.
Stability	
Store Temp	Ambient
Ship Temp	
Product Description	A major metabolites of benzene. It causes DNA damage and apoptosis through H ₂ O ₂ generation in cells.
References	Hiraku Y., Kawanishi S. <i>Can. Res.</i> 56 :5172-8 (1996).



Pricing

100 g	\$24.30
500 g	\$68.90

Bulk quantities available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
St. Paul, MN 55114
USA

Phone: 888-LKT-LABS (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.