




LKT Laboratories, Inc
2233 University Avenue West
St. Paul, Minnesota 55114 USA
Ph: 651.644.8424 Fax: 651.644.8357
www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	90357-06-5
Product ID	B3209
Product Name	Bicalutamide
Chemical Name	N-(4-cyano-3-(trifluoromethyl)phenyl)-3-((4-fluorophenyl)sulfonyl)-2-hydroxy-2-methyl-, (+)-Propanamide
Synonym	Casodex
Formula	C ₁₈ H ₁₄ F ₄ N ₂ O ₄ S
Formula Wt.	430.37
Melting Point	191°C-193°C
Boiling Point	
Purity	≥98%
Solubility	Soluble in acetone or THF.
Stability	
Store Temp	
Ship Temp	
Product Description	A nonsteroidal antiandrogen, which has been shown effective in the treatment of prostate cancer. It has also been shown to induce cell death by a pathway that is independent of changes in mitochondrial membrane potential and Bcl-2 actions.
References	Maucher A, von Angerer E. J Canc Res Clin Oncol. 119:669-74 (1993). Lee EC, Zhan P, Schallhom R et al. Cell Death & Differentiation. 10:761-71 (2003).



Pricing

100 mg	\$50.20
250 mg	\$90.90
1 g	\$273.50

Bulk quantities available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
St. Paul, MN 55114
USA

Phone: 888-LKT-LABS (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.