




LKT Laboratories, Inc
2233 University Avenue West
St. Paul, Minnesota 55114 USA
Ph: 651.644.8424 Fax: 651.644.8357
www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	113299-40-4
Product ID	B3210
Product Name	R-Bicalutamide
Chemical Name	R-N-(4-cyano-3-(trifluoromethyl)phenyl)-3-((4-fluorophenyl)sulfonyl)-2-hydroxy-2-methyl-, (+)-Propanamide
Synonym	
Formula	C ₁₈ H ₁₄ F ₄ N ₂ O ₄ S
Formula Wt.	430.37
Melting Point	178-179°C
Boiling Point	
Purity	
Solubility	water: 5mg/L
Stability	
Store Temp	
Ship Temp	
Product Description	[a] _D ²⁰ = -83.~ - 85.0 (C=1, MeOH)
References	



Pricing

100 mg	\$100.40
250 mg	\$181.90
1 g	\$547.00

Bulk quantities available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
St. Paul, MN 55114
USA

Phone: 888-LKT-LABS (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.