




LKT Laboratories, Inc
2233 University Avenue West
St. Paul, Minnesota 55114 USA
Ph: 651.644.8424 Fax: 651.644.8357
www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	
Product ID	B3280
Product Name	Bis(salicyl) fumarate
Chemical Name	
Synonym	
Formula	C ₁₈ H ₁₂ O ₈
Formula Wt.	356.29
Melting Point	178°C-180°C
Boiling Point	
Purity	≥98%
Solubility	
Stability	
Store Temp	-20 deg C
Ship Temp	Ambient
Product Description	A potent acylating agents of intracellular hemoglobin. It cross-links beta chains of hemoglobin.
References	Zaugg, R.H.; Walder, J.A.; Walder, R.Y.; Steele, J.M.; Klotz, I.M. <i>J. Biol. Chem.</i> 255 :2816-21 (1980).



Pricing

100 mg	\$16.50
500 mg	\$57.50
1 g	\$82.00
5 g	\$164.00

Bulk quantities available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
St. Paul, MN 55114
USA

Phone: 888-LKT-LABS (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.