




LKT Laboratories, Inc
2233 University Avenue West
St. Paul, Minnesota 55114 USA
Ph: 651.644.8424 Fax: 651.644.8357
www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	33570-04-6
Product ID	B3345
Product Name	(-)-Bilobalide
Chemical Name	(5aR-(3aS,5aα,8b,8aS,9a,10aα))-9-(1,1-dimethylethyl)-10,10a-dihydro-8,9-dihydroxy-4H,5aH,9H-furo[2,3-b]furo[3',2':2,3]cyclopenta[1,2-c]furan-2,4,7(3H,8H)-trione
Synonym	Bilobalide
Formula	C ₁₅ H ₁₈ O ₈
Formula Wt.	326.30
Melting Point	
Boiling Point	
Purity	≥98%
Solubility	Soluble in ethanol or
Stability	
Store Temp	-20 deg C
Ship Temp	
Product Description	Bilobalide is a naturally occurring diterpene from <i>Ginkgo biloba</i> . It has anti-ischemic and anti-convulsant properties.
References	Janssens, D., Remacle, J., Drieu, K. Michiels, C. <i>Biochem. Pharmacol.</i> 58 :109-19 (1999). Sasaki, K., Hatta, S., Haga, M., Ohshika, H., <i>Eur. J. Pharmacol.</i> 367 :165-73 (1999).



Pricing

5 mg	\$112.10
10 mg	\$172.10
25 mg	\$344.40

Bulk quantities available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
St. Paul, MN 55114
USA

Phone: 888-LKT-LABS (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.