




LKT Laboratories, Inc
2233 University Avenue West
St. Paul, Minnesota 55114 USA
Ph: 651.644.8424 Fax: 651.644.8357
www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	13687-09-7
Product ID	B3373
Product Name	Bis(aziridinyl)methylamino phosphine sulfide / Bis(aziridinyl)methylamino phosphine sulfide
Chemical Name	Bis(aziridinyl)methylamino phosphine sulfide
Synonym	Bisazir
Formula	C ₃ H ₁₂ N ₃ SP
Formula Wt.	177.21
Melting Point	69°C-72°C
Boiling Point	
Purity	≥98%
Solubility	Soluble in DMSO.
Stability	
Store Temp	-20 deg C
Ship Temp	Ambient
Product Description	Effects male sterility in Anopheles albimanus mosquitoes.
References	



Pricing

500 mg	\$59.60
1 g	\$100.70
5 g	\$400.30

Bulk quantities available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
St. Paul, MN 55114
USA

Phone: 888-LKT-LABS (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.