




LKT Laboratories, Inc

2233 University Avenue West
St. Paul, Minnesota 55114 USA
Ph: 651.644.8424 Fax: 651.644.8357
www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	603-50-9
Product ID	B3374
Product Name	Bisacodyl
Chemical Name	4,4'-(2-Pyridylmethylene)-bisphenol diacetate
Synonym	Bicol; Broxalax; Contalax; DAMP; Dulrolax; Endokolat; Laxadin; Laxorex; Perilax; Teletin; Ulcolax
Formula	C ₂₂ H ₁₉ NO ₄
Formula Wt.	361.39
Melting Point	138°C
Boiling Point	
Purity	≥98%
Solubility	Soluble in acids, acetone,
Stability	
Store Temp	+4 deg C
Ship Temp	
Product Description	A diphenolic laxative. It is suggested that intestinal inhibition of (Na+K) ATPase activity, increase of mucosal PGE2 content, and possibly also stimulation of adenylyl cyclase activity might contribute to the net water accumulation induced by bisacodyl.
References	Farack UM, Nell G. <i>Digestion</i> . 30 :191-4 (1984). Rachmilewitz D, Karmeli F, Okon E. <i>Dig Dis Sci</i> . 25 :602-8 (1980).



Pricing

10 g	\$66.30
25 g	\$71.80

Bulk quantities available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
St. Paul, MN 55114
USA

Phone: 888-LKT-LABS (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.