




LKT Laboratories, Inc
2233 University Avenue West
St. Paul, Minnesota 55114 USA
Ph: 651.644.8424 Fax: 651.644.8357
www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	22150-76-1
Product ID	B3458
Product Name	Biopterin
Chemical Name	2-Amino-6-[(1 <i>R</i> ,2 <i>S</i>)-1,2-dihydroxypropyl]-4(1 <i>H</i>)- pteridinone
Synonym	Pterin HB ₂
Formula	C ₉ H ₁₁ N ₅ O ₃
Formula Wt.	237.22
Melting Point	250°C-280°C (chars)
Boiling Point	
Purity	≥98%
Solubility	Soluble in alcohol, ether,
Stability	
Store Temp	Ambient
Ship Temp	Ambient
Product Description	An obligate cofactor of inducible nitric oxide synthase.
References	Weinberg, J.B., Misukonis, M.A., Shami, P.J. <i>et. al. Blood.</i> 86 :1184-95 (1995).



Pricing

5 mg	\$100.80
10 mg	\$168.00
25 mg	\$358.40

Bulk quantities available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
St. Paul, MN 55114
USA

Phone: 888-LKT-LABS (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.