




**LKT Laboratories, Inc**  
2233 University Avenue West  
St. Paul, Minnesota 55114 USA  
Ph: 651.644.8424 Fax: 651.644.8357  
www.lktlabs.com info@lktlabs.com

# Product Information

<b>CAS No.</b>	55658-47-4
<b>Product ID</b>	B4517
<b>Product Name</b>	Bleomycin A5 hydrochloride
<b>Chemical Name</b>	
<b>Synonym</b>	Pingyangmycin hydrochloride
<b>Formula</b>	
<b>Formula Wt.</b>	1440.6
<b>Melting Point</b>	
<b>Boiling Point</b>	
<b>Purity</b>	≥90%
<b>Solubility</b>	
<b>Stability</b>	
<b>Store Temp</b>	+4 deg C
<b>Ship Temp</b>	
<b>Product Description</b>	Bleomycin A5 hydrochloride is an antitumor antibiotic which induces apoptosis.
<b>References</b>	Tai, K.W., Chou, M.Y., Hu, C.C., Yang, J.J., Chang, Y.C. <i>Oral Oncol.</i> <b>36</b> :242-7 (2000). Li, X.T. <i>Chung Kuo I Hsueh Ko Hsueh Yuan Hsueh Pao.</i> <b>12</b> :182-6 (1990).



## Pricing

5 mg	\$261.70
10 mg	\$438.90
25 mg	\$903.70

**Bulk quantities available upon request**

## LKT Laboratories, Inc.

Address: 2233 University Ave. W.  
St. Paul, MN 55114  
USA

Phone: 888-LKT-LABS (558-5227)  
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

**Caution:** This product is intended for laboratory and research use only. It is not for human or drug use.