



LKT Laboratories, Inc
 2233 University Avenue West
 St. Paul, Minnesota 55114 USA
 Ph: 651.644.8424 Fax: 651.644.8357
 www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	31362-50-2
Product ID	B5648
Product Name	Bombesin
Chemical Name	
Synonym	2-L-glutamine-6-L-asparagine-Alytesin
Formula	C ₇₁ H ₁₁₁ N ₂₄ O ₁₈ S
Formula Wt.	1619.86
Melting Point	
Boiling Point	
Purity	≥97%
Solubility	Soluble in 5% acetic acid.
Stability	Hygroscopic.
Store Temp	+4 deg C
Ship Temp	



Product Description A gut tetradecapeptide with the ability to stimulate the release of numerous hormones. It inhibits growth of pancreatic ductal adenocarcinoma (H2T) in nude mice.

References Chen LW, Hsu CM, Huang JK et al. J Formos Med Assoc. 99:491-8 (2000).
 Farre A, Ishizuka J, Gomez G et al. Pancreas. 9:652-6 (1994).

Pricing	
1 mg	\$28.70
2 mg	\$48.40
5 mg	\$86.00

Bulk quantities available upon request

LKT Laboratories, Inc.
 Address: 2233 University Ave. W.
 St. Paul, MN 55114
 USA
 Phone: 888-LKT-LABS (558-5227)
 651-644-8424
 Fax: 651-644-8357
 Email: info@lktlabs.com
 Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.