




LKT Laboratories, Inc
2233 University Avenue West
St. Paul, Minnesota 55114 USA
Ph: 651.644.8424 Fax: 651.644.8357
www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	11076-19-0
Product ID	B5753
Product Name	Bongkreikic acid
Chemical Name	3-Carboxymethyl-17-methoxy-6,18,21-trimethyldocosa-2,4,8,12,14,18,20-heptaenedioic acid
Synonym	
Formula	C ₂₈ H ₃₈ O ₇
Formula Wt.	486.60
Melting Point	50°C-60°C
Boiling Point	
Purity	≥98%
Solubility	Soluble as 2N in NH ₄ OH.
Stability	Protect from light.
Store Temp	-20 deg C
Ship Temp	
Product Description	A mitochondrial permeability transition pore blocker. It has been shown to inhibit apoptosis by preventing PARP cleavage and DEVDase activity.
References	Yoon HS, Moon SC, Kim ND et al. Biochem Bioph Res Co. 276:151-6 (2000). Zamora M, Granell M, Mampel T et al. FEBS Letters. 563:155-60 (2004).



Pricing

100 µg	\$307.30

Bulk quantities available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
St. Paul, MN 55114
USA

Phone: 888-LKT-LABS (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.