



Product Information

CAS Registry No. 80451-05-4

Product Name: Cecropin B

Catalog No. C1609

Chemical Name:

Synonym(s):

Mol. formula: $C_{176}H_{302}N_{52}O_{41}S$

Formula wt: 3834.76

M.p.

B.p.

Purity:

Solubility:

Stability (in original shipping container):

Storage temp:

Product Description:

Antibacterial peptide isolated from pig intestine and moths.

H-Lys-Trp-Lys-Val-Phe-Lys-Lys-Ile-Glu-Lys-Met-Gly-Arg-Asn-Ile-Arg-Asn-Gly-Ile-Val-Lys-Ala-Gly-Pro-Ala-Ile-Ala-Val-Leu-Gly-Glu-Ala-Lys-Ala-Leu-NH₂

References:

Lee, J.-Y., et al. Proc. Nat. Acad. Sci. USA 86:9159-9162 (1989).

PRICING

1 mg	\$185.60
2 mg	\$315.20
5 mg	\$556.80

Bulk quantity available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.