




LKT Laboratories, Inc
2233 University Avenue West
St. Paul, Minnesota 55114 USA
Ph: 651.644.8424 Fax: 651.644.8357
www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	50-41-9
Product ID	C4559
Product Name	Clomiphene Citrate
Chemical Name	2-[4-(2-Chloro-1,2-diphenylethenyl)phenoxy]-N,N-diethylethanamine
Synonym	Clomid; Clomphid; Clostilbegyt; Dyneric; Pergotime; Serophene
Formula	C ₂₆ H ₂₈ ClNO.C ₆ H ₈ O ₇
Formula Wt.	598.09
Melting Point	116.5°C-118°C
Boiling Point	
Purity	≥98%
Solubility	Slightly soluble in water and
Stability	
Store Temp	+4 deg C
Ship Temp	
Product Description	A selective oestrogen receptor modulator. Used as an antileukemic drug which kill cells by apoptosis mediated by oxidative stress and activation of PKC.
References	Haskell SG. <i>South Med J.</i> 96 :469-76 (2003). Hayon T, Dvilansky A, Oriev L, Nathan I. <i>Anticancer Res.</i> 19 :2089-93 (1999).



Pricing

	1 g	\$27.70
	5 g	\$75.90
	10 g	\$124.20

Bulk quantities available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
St. Paul, MN 55114
USA

Phone: 888-LKT-LABS (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.