



Product Information

CAS Registry No. 66640-93-5

Product Name: Difluoromethylornithine

Catalog No. D3221

Chemical Name: 2-(Difluoromethyl)-DL-ornithine

Synonym(s): Eflornithine; DFMO

Mol. formula: $C_6H_{12}F_2N_2O_2 \cdot HCl \cdot H_2O$

Formula wt: 236.65

M.p.

B.p.

Purity: $\geq 98\%$


Solubility: Soluble in water. Insoluble in acetone or chloroform.

Stability (in original shipping container):

Storage temp: -20 deg C

Product Description:

An irreversible inhibitor of ornithine decarboxylase.



References:

Verma, A.K. *Prev.Med.* **18**:646-652 (1989).

PRICING

10 mg \$46.10

25 mg \$109.70

*Bulk quantity available
upon request*

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.