



Product Information

CAS Registry No. 587-63-3

Product Name: 7,8-Dihydrokawain

Catalog No. D3229

Chemical Name:

Synonym(s): DHK; Dihydrokawain; Marindinin

Mol. formula: $C_{14}H_{16}O_3$

Formula wt: 232.28

M.p.


B.p.

Purity: $\geq 98\%$

Solubility: Soluble in DMSO and methanol.

Stability (in original shipping container):

Storage temp: -20 deg C



Product Description:

One of the six main active ingredients (a-pyrone) found in kava-kava, exhibits neuroprotective and antinociceptive activity.

References:

Backhauss, C., Krieglstein, J. *Eur J Pharmacol* **215**:265-9 (1992).

Jamieson, D.D., Duffield P.H. *Clin Exp Pharmacol Physiol.*

17:495-507 (1990).

PRICING

5 mg	\$99.50
10 mg	\$153.70

*Bulk quantity available
upon request*

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.