



Product Information

CAS Registry No. 22494-42-4

Product Name: **Diflunisal**

Catalog No. **D3322**

Chemical Name: 2',4'-Difluoro-4-hydroxy-[1,1'-biphenyl]-3-carboxylic acid

Synonym(s): Adomal; Difludol; Dolisal; Dolobid; Dolobis; Flovacil; Fluniget; Fluodonil

Mol. formula: $C_{13}H_8F_2O_3$

Formula wt: 250.20

M.p. 210°C-211°C

B.p.

Purity: $\geq 98\%$


Solubility: Sparingly soluble in water.

Stability (in original shipping container):

Storage temp:

Product Description:

A non-steroidal anti-inflammatory analgesic that has unusually long duration of action. Recently, it was found to have anti-proliferative activity against colon adenocarcinoma cells.



References:

Forbes, JA., Beaver, WT., White, EH., *et al.* *JAMA*. **248**:2139-2142 (1982).

Cannell, GR., Vesey, DA., Dickinson, RG. *Life Sci.* **70**:37-48 (2001).

PRICING

5 g	\$34.00
10 g	\$56.90
50 g	\$203.30

Bulk quantity available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.