



Product Information

CAS Registry No. 24512-63-8

Product Name: Geniposide

Catalog No. G1650

Chemical Name:

Synonym(s):

Mol. formula: $C_{17}H_{24}O_{10}$

Formula wt: 388.37

M.p. 163°C-164°C

B.p.

Purity: $\geq 98\%$


Solubility: Soluble in water, methanol and ethanol.

Stability (in original shipping container):

Storage temp:

Product Description:

Natural product isolated from the fruit of Gerdenia. It was found to induce increased activity of phase II detoxifying enzymes, inhibit tumor promotion, and induce apoptosis in rat C6 glioma cells.



References:

Wang, C.J., Wang, S.W., Lin, J.K. *Cancer Letters* **60**:95-102 (1991).
Lee, M.J., Hsu, J.D., Wang, C.J. *Anticancer Res.* **15**:411-416 (1995).
Chang, Y.C., Tseng, T.H., Lee, M.J. *Chemico-Biol Interactions* **141**: 243-257 (2002).

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.

PRICING

10 mg	\$78.70
25 mg	\$169.50
100 mg	\$542.10

Bulk quantity available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com