



# Product Information

**CAS Registry No.** 89468-62-2

**Product Name:** **Helodormin**

**Catalog No.** **H1644**

**Chemical Name:** 5-L-isoleucine-8-L-glutamine-9-L-glutamine-24-L-alanine-30-L-arginine-34-de-L-arginine-36-L-prolinamide-37-de-L-serine-38-de-L-serine- Helospectin I

**Synonym(s):** Exendin 2; Helodermin; Heliodermin.

**Mol. formula:**  $C_{176}H_{285}N_{47}O_{49}$

**Formula wt:** 3843.47

**M.p.**

**B.p.**

**Purity:**  $\geq 98\%$

**Solubility:**

**Stability** (in original shipping container):

**Storage temp:** -20 deg C

**Product Description:**

His-Ser-Asp-Ala-Ile-Phe-Thr-Gln-Gln-Tyr-Ser-Lys-Leu-Leu-Ala-Lys-Leu-Ala-Leu-Gln-Lys-Tyr-Leu-Ala-Ser-Ile-Leu-Gly-Ser-Arg-Thr-Ser-Pro-Pro-Pro-NH<sub>2</sub>

**References:**

**Caution:** This product is intended for laboratory and research use only. It is not for human or drug use.

## PRICING

1 mg \$368.00

*Bulk quantity available upon request*

## LKT Laboratories, Inc.

**Address:** 2233 University Ave. W  
St. Paul, MN 55114  
USA

**Phone:** 888-LKT-Labs (558-5227)  
651 644 8424

**Fax:** 651-644-8357

**Email:** info@lktlabs.com

**Web:** www.lktlabs.com