



# Product Information

**CAS Registry No.**

**Product Name:** Hemorphin-7

**Catalog No.** H1648

**Chemical Name:**

**Synonym(s):**

**Mol. formula:**  $C_{49}H_{64}N_{12}O_{11}$

**Formula wt:** 997.12

**M.p.**

**B.p.**

**Purity:**  $\geq 98\%$

**Solubility:**

**Stability** (in original shipping container):

**Storage temp:** -20 deg C

**Product Description:**

Tyr-Pro-Trp-Thr-Gln-Arg-Phe

**References:**

## PRICING

1 mg	\$64.00
2 mg	\$108.80
5 mg	\$192.00

*Bulk quantity available  
upon request*

## LKT Laboratories, Inc.

**Address:** 2233 University Ave. W  
St. Paul, MN 55114  
USA

**Phone:** 888-LKT-Labs (558-5227)  
651 644 8424

**Fax:** 651-644-8357

**Email:** info@lktlabs.com

**Web:** www.lktlabs.com

**Caution:** This product is intended for laboratory and research use only. It is not for human or drug use.