



Product Information

CAS Registry No. 520-33-2

Product Name: Hesperetin

Catalog No. H1672

Chemical Name: (2S)-2,3-Dihydro-5,7-dihydroxy-2-(3-hydroxy-4-methoxyphenyl)-4H-1-benzopyran-4-one

Synonym(s): 3',5,7-trihydroxy-4'-methoxyflavanone
; Cyanidanon 4'-methyl ether 1626

Mol. formula: C₁₆H₁₄O₆

Formula wt: 302.28

M.p. 226°C-228°C

B.p.

Purity: ≥95%


Solubility: Soluble in ethanol. Slightly soluble in water.

Stability (in original shipping container):

Storage temp: Ambient

Product Description:

The aglycone of Hesperidin.



References:

PRICING

1 g	\$18.40
5 g	\$57.40
10 g	\$91.60

*Bulk quantity available
upon request*

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.