



Product Information

CAS Registry No. 52806-53-8

Product Name: 2-Hydroxy-flutamide

Catalog No. H9718

Chemical Name:

Synonym(s): Hydroxyniphtholide

Mol. formula: C₁₁H₁₁F₃N₂O₄

Formula wt: 292.21

M.p.

B.p.

Purity: ≥98%


Solubility:

Stability (in original shipping container):

Storage temp: -20 deg C

Product Description:

Metabolite of flutamide found to inhibit the proliferation of estradiol-induced growth of MCF-7 breast cancer cells.



References:

Shet, M.S., McPhaul, M., Fisher, CW., *et al Drug Metab Dispos.* **25**:1298-303 (1997).

Di Monaco, M., Brignardello, E., Leonardi L., *et al J Cancer Res Clin Oncol* . **121**:710-4 (1995).

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.

PRICING

10 mg	\$95.10
25 mg	\$183.10
100 mg	\$658.70

Bulk quantity available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com