



Product Information

CAS Registry No. 3155-48-4

Product Name: **Kawain**

Catalog No. **K0088**

Chemical Name: (6*R*)-5,6-Dihydro-4-methoxy-6-[(1*E*)-2-phenyl-ethenyl]-2*H*-pyran-2-one

Synonym(s): Kavain; Gonosan

Mol. formula: C₁₄H₁₄O₃

Formula wt: 230.26

M.p. 105°C-106°C

B.p. 195°C-197°C

Purity: ≥98%


Solubility: Soluble in acetone, ether or methanol. Slightly soluble in hexane. Practically insoluble in water.

Stability (in original shipping container):

Storage temp: -20 deg C

Product Description:

A major component of kava kava extract, has antinociceptive activity.



References:

Jamieson, DD., Duffield PH. *Clin Exp Pharm Physiol.* **16**:496-507 (1990).

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.

PRICING

5 mg	\$99.50
10 mg	\$153.70

Bulk quantity available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com