



Product Information

CAS Registry No.

Product Name: Kemptide

Catalog No. K1650

Chemical Name:

Synonym(s):

Mol. formula: $C_{32}H_{61}N_{13}O_9$

Formula wt: 771.92

M.p.

B.p.

Purity:

Solubility:

Stability (in original shipping container):

Storage temp:

Product Description:

A synthetic heptapeptide substrate for kinases that maintains cell membrane intactness when added to cultured cells.

H-Leu-Arg-Arg-Ala-Ser-Leu-Gly-OH

References:

Foxwell BM, Band HA, Long J et al. Br J Cancer. 57:489-493 (1988).

Kubler D, Pyerin W, Bill O et al. J Biol Chem. 264:14549-14555 (1989).

PRICING

1 mg	\$32.00
2 mg	\$54.40
5 mg	\$96.00

Bulk quantity available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.