



Product Information

CAS Registry No.

Product Name: Kinetensin

Catalog No. K3352

Chemical Name:

Synonym(s):

Mol. formula: $C_{56}H_{85}N_{17}O_{11}$

Formula wt: 1172.4

M.p.

B.p.

Purity:

Solubility:

Stability (in original shipping container):

Storage temp:

Product Description:

Increases vascular permeability and induces histamine release in rat peritoneal mast cells.

H-Ile-Ala-Arg-Arg-His-Pro-Tyr-
Phe-Leu-OH

References:

Sydbom A, Ware J, Mogard MH. Agents Actions. 27:68-71 (1989).

PRICING

1 mg	\$32.00
2 mg	\$54.40
5 mg	\$96.00

*Bulk quantity available
upon request*

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.