



# Product Information

**CAS Registry No.**

**Product Name:** M40

**Catalog No.** M0040

**Chemical Name:**

**Synonym(s):**

**Mol. formula:** C<sub>94</sub>H<sub>145</sub>N<sub>23</sub>O<sub>24</sub>

**Formula wt:** 1981.34

**M.p.**

**B.p.**

**Purity:**

**Solubility:**

**Stability** (in original shipping container):

**Storage temp:**

**Product Description:**

A galanin receptor antagonist.

H-Gly-Trp-Thr-Leu-Asn-Ser-Ala-  
Gly-Tyr-Leu-Leu-Gly-Pro-Pro-  
Pro-Ala-Leu-Ala-Leu-Ala-NH<sub>2</sub>

**References:**

Crawley JN, Robinson JK, Langel U, Bartfai T. Brain Research; 600(2): 268-272 (1993).

## PRICING

0.5 mg	\$121.60
1 mg	\$206.40
2.5 mg	\$364.80

*Bulk quantity available  
upon request*

## LKT Laboratories, Inc.

**Address:** 2233 University Ave. W  
St. Paul, MN 55114  
USA

**Phone:** 888-LKT-Labs (558-5227)  
651 644 8424

**Fax:** 651-644-8357

**Email:** info@lktlabs.com

**Web:** www.lktlabs.com

**Caution:** This product is intended for laboratory and research use only. It is not for human or drug use.