



Product Information

CAS Registry No. 1006-94-6

Product Name: 5-Methoxyindole

Catalog No. M1680

Chemical Name: 5-Methoxy-1H-indole

Synonym(s): Femedol

Mol. formula: C₉H₉NO

Formula wt: 147.17

M.p. 52°C-55°C

B.p. 176°C-178°C 17mm Hg (lit.)


Purity: ≥98%

Solubility: Insoluble in water.

Stability (in original shipping container):

Storage temp: +4 deg C

Product Description:



References:

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.

PRICING

1 g	\$52.50
5 g	\$207.50

*Bulk quantity available
upon request*

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com