



Product Information

CAS Registry No. 34807-41-5

Product Name: Mezerein

Catalog No. M1699

Chemical Name: 12β-[(E,E)-5-Phenyl-2,4-pentadienyloxy]daphnetoxin

Synonym(s): Meserein

Mol. formula: C₃₈H₃₈O₁₀

Formula wt: 654.70

M.p. 256°C-269°C

B.p.

Purity: ≥98%


Solubility: Soluble in DMSO or ethanol.

Stability (in original shipping container):

Storage temp: -20 deg C

Product Description:

A protein kinase C activator, second stage tumor promoter.



References:

Slaga, T.J., *et al.* *Proc. Nat. Acad. Sci. USA* **77**:3659-3663 (1980).

PRICING

1 mg	\$97.60
5 mg	\$338.80
10 mg	\$609.90

*Bulk quantity available
upon request*

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.