



Product Information

CAS Registry No. 22916-47-8

Product Name: Miconazole

Catalog No. M3309

Chemical Name: 1-[2-(2,4-Dichlorophenyl)-2-[(2,4-dichlorophenyl)-methoxy]ethyl]-1H-imidazole

Synonym(s): Monistat IV

Mol. formula: C₁₈H₁₄Cl₄N₂O

Formula wt: 416.14

M.p. 84.3°C


B.p.

Purity: ≥98%

Solubility: Practically insoluble in water.

Stability (in original shipping container):

Storage temp: Ambient



Product Description:

An antifungal agent with good aromatase inhibitory activity, found to be an effective chemopreventive in a rat mammary cancer model.

References:

Lubert, R.A., Steele, V.E., Etyo, I., *et al.* *Int. J. Cancer.* **72**:95-101 (1997).

PRICING

1 g	\$23.80
5 g	\$82.70
25 g	\$265.80

Bulk quantity available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.