



Product Information

CAS Registry No. 22204-53-1

Product Name: D-Naproxen

Catalog No. N0061

Chemical Name: S(+)-6-Methoxy- α methyl-2-naphthaleneacetic acid

Synonym(s): Nonyl; Equiproxen; Napren; Pranoxen

Mol. formula: $C_{14}H_{14}O_3$

Formula wt: 230.26

M.p. 155°C-158°C

B.p.

Purity: $\geq 98\%$


Solubility: Soluble in ethanol and acetone. Insoluble in water.

Stability (in original shipping container):

Storage temp: Ambient

Product Description:

A non-steroidal anti-inflammatory agent found to significantly inhibit NNK activation in lung tissues.



References:

Bouchard, L., Castonguay, A. *Drug Metab. Disp.* **21**:293-298 (1993).

PRICING

5 g	\$35.20
25 g	\$122.90
50 g	\$224.50

*Bulk quantity available
upon request*

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.