



Product Information

CAS Registry No. 22204-53-1

Product Name: D,L-Naproxen

Catalog No. N0062

Chemical Name: (+)-(S)-6-Methoxy- α -methyl-2-naphthaleneacetic acid

Synonym(s): Artagen; Danaprox; Nafasol; Sutolin

Mol. formula: C₁₄H₁₄O₃

Formula wt: 230.26

M.p. 152°C-154°C


B.p.

Purity: ≥96%

Solubility: Insoluble in water. Soluble in acetone or ethanol.

Stability (in original shipping container):

Storage temp: Ambient



Product Description:

A non-steroidal anti-inflammatory agent found to significantly inhibit NNK activation in lung tissues.

References:

Bouchard, L., Castonguay, A. *Drug Metab.Disp.* **21**:293-298 (1993).

PRICING

10 g	\$36.10
25 g	\$80.90
100 g	\$256.00

*Bulk quantity available
upon request*

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.