



Product Information

CAS Registry No. [66085-59-4]

Product Name: Nimodipine

Catalog No. N3448

Chemical Name: 1,4-Dihydro-2,6-dimethyl-4-(3-nitrophenyl)-3,5-pyridinedicarboxylic acid 2-methoxyethyl 1-methyl-ethyl ester

Synonym(s): Periplum; Nimotop

Mol. formula: $C_{21}H_{26}N_2O_7$

Formula wt: 418.44

M.p. 125°C

B.p.

Purity: ≥98%


Solubility: Soluble in methanol or DMSO. Insoluble in water.

Stability (in original shipping container):

Storage temp:

Product Description:

A calcium channel antagonist.



References:

Gupta MC, Garg SK, Das BP, Bhargava VK. Indian J Physiol Pharmacol. 47:347-51 (2003).

PRICING

| | |
|--------|----------|
| 500 mg | \$49.30 |
| 1 g | \$86.30 |
| 5 g | \$308.00 |

*Bulk quantity available
upon request*

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.