



Product Information

CAS Registry No. 51803-78-2

Product Name: Nimesulide

Catalog No. N3450

Chemical Name: N-(4-Nitro-2-phenoxyphenyl)methanesulfonamide

Synonym(s): Aulin; Fansidol; Flogovital; Mesulid; Nide; Nidol; Nisulid

Mol. formula: C₁₃H₁₂N₂O₅S

Formula wt: 308.31

M.p. 143°C-144.5°C

B.p.

Purity: ≥98%


Solubility: Soluble in chloroform.

Stability (in original shipping container):

Storage temp: Ambient

Product Description:

Nimesulide inhibits prostaglandin synthetase and exhibits anti-inflammatory activity.



References:

Vigdahl, R.; Tukey, R. *Biochem. Pharmacol.* **26**:307 (1977).
Swingle, K.; More, G. *Drug Exp. Clin. Res.* **10**:587 (1984).

PRICING

5 g	\$38.50
10 g	\$69.30
25 g	\$138.40

*Bulk quantity available
upon request*

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.