



Product Information

CAS Registry No. 170713754

Product Name: Nociceptin

Catalog No. N5210

Chemical Name: Orphanin FQ

Synonym(s):

Mol. formula: $C_{79}H_{129}N_{27}O_{22}$

Formula wt: 1809.1

M.p.

B.p.

Purity: $\geq 98\%$

Solubility: Soluble in water.

Stability (in original shipping container):

Storage temp: -20 deg C

Phe-Gly-Gly-Phe-Thr-Gly-Ala-Arg-Lys-Ser-Ala-Arg-Lys-Leu-Ala-Asn-Gln

Product Description:

The endogenous ligand for the opioid-receptor-like receptor 1. It mediates essential functions in the central and peripheral nervous systems.

References:

Li N, Wei SY, Yu LC et al. Brain Res. 1025:67-74 (2004).

PRICING

1 mg \$161.30

Bulk quantity available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.