



Product Information

CAS Registry No. 31430-18-9

Product Name: Nocodazole

Catalog No. N5409

Chemical Name: [5-(2-Thienylcarbonyl)-1H-benzimidazol-2-yl]- carbamic acid methyl ester

Synonym(s): Oncodazole; R 17934

Mol. formula: $C_{14}H_{11}N_3O_3S$

Formula wt: 301.32

M.p. 300°C (dec)

B.p.

Purity: $\geq 98\%$


Solubility: Soluble in DMSO. Insoluble in water.

Stability (in original shipping container):

Storage temp: -20 deg C

Product Description:

Nocodazole is a microtubule-active drug.



References:

Karbowski M, Spodnik JH, Teranishi Ma, *et al. J. Cell Sci.* **114**:281-291 (2001).

PRICING

10 mg	\$95.00
50 mg	\$350.00

*Bulk quantity available
upon request*

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.