



Product Information

CAS Registry No. 6833-84-7

Product Name: Nonactin

Catalog No. N5652

Chemical Name:

Synonym(s):

Mol. formula: $C_{40}H_{64}O_{12}$

Formula wt: 736.934

M.p.

B.p.

Purity:


Solubility:

Stability (in original shipping container):

Storage temp:

Product Description:

A selective ionophore used in the determination of ammonium ions.



References:

Karakus E, Pekyardimci S, Kilic E. Artif Cells Blood Substit Immobil Biotechnol. 34:523-34 (2006).

PRICING

| | |
|-------|---------|
| 1 mg | \$20.00 |
| 5 mg | \$40.00 |
| 10 mg | \$72.00 |

*Bulk quantity available
upon request*

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.