



Product Information

CAS Registry No. 70458-96-7

Product Name: Norfloxacin

Catalog No. N5768

Chemical Name: 1-Ethyl-6-fluoro-1,4-dihydro-4-oxo-7-(1-piperazinyl)-3-quinolinecarboxylic acid

Synonym(s): Baccidal; Sebercim

Mol. formula: $C_{16}H_{18}FN_3O_3$

Formula wt: 319.33

M.p. 227-228°C

B.p.

Purity:


Solubility: Soly at 25°C (mg/ml): water 0.28; methanol 0.98; ethanol 1.9.

Stability (in original shipping container):

Storage temp: Ambient

Product Description:

A fluorinated antibacterial quinolone.



References:

PRICING

10 g	\$56.20
50 g	\$199.70

*Bulk quantity available
upon request*

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.