



Product Information

CAS Registry No. 57186-25-1

Product Name: Paxilline

Catalog No. P0392

Chemical Name:

Synonym(s):

Mol. formula: $C_{27}H_{33}NO_4$

Formula wt: 435.56

M.p.

B.p.

Purity:


Solubility:

Stability (in original shipping container):

Storage temp:

Product Description:

A tremorgenic indole alkaloid that selectively blocks high-conductance Ca^{2+} -activated K^{+} -channels.



References:

Sanchez M et al Neuropharmacology 35:963 (1996).

Young C et al Curr. Genet. 33:368 (1998).

PRICING

5 mg	\$118.00
10 mg	\$216.00

*Bulk quantity available
upon request*

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.