




**LKT Laboratories, Inc**  
2233 University Avenue West  
St. Paul, Minnesota 55114 USA  
Ph: 651.644.8424 Fax: 651.644.8357  
www.lktlabs.com info@lktlabs.com

# Product Information

<b>CAS No.</b>	87333-19-5
<b>Product ID</b>	R0249
<b>Product Name</b>	Ramipril
<b>Chemical Name</b>	[2S,3aS,6aS]-1-[(2S)-2-[[[(1S)-1-(Ethoxycarbonyl)-3-phenylpropyl]amino]-1-oxopropyl]octahydrocyclopenta[b]pyrrole-2-carboxylic acid
<b>Synonym</b>	Tritace; Ramace; Unipril; Vesdil
<b>Formula</b>	C <sub>23</sub> H <sub>32</sub> N <sub>2</sub> O <sub>5</sub>
<b>Formula Wt.</b>	416.51
<b>Melting Point</b>	109°C
<b>Boiling Point</b>	
<b>Purity</b>	≥98%
<b>Solubility</b>	Soluble in DMSO
<b>Stability</b>	
<b>Store Temp</b>	
<b>Ship Temp</b>	
<b>Product Description</b>	An angiotensin-converting enzyme (ACE) inhibitor. Ramipril is a prodrug of the active metabolite ramiprilat which is formed upon absorption.
<b>References</b>	Warner GT, Perry CM. Am J Cardiovasc Drugs. 3:113-6 (2003).



## Pricing

500 mg	\$69.00
1 g	\$107.60
5 g	\$414.00

**Bulk quantities available upon request**

## LKT Laboratories, Inc.

Address: 2233 University Ave. W.  
St. Paul, MN 55114  
USA

Phone: 888-LKT-LABS (558-5227)  
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

**Caution:** This product is intended for laboratory and research use only. It is not for human or drug use.